INFORMATION LEAFLET
NOTICE EXPLICATIVE
PACKUNGSBEILAGE
INDLÆGSSEDDEL
BIJSLUITER
FOLLETO ILLUSTRATIVO
FOGLIO ILLUSTRATIVO
INFORMATIONSBLAD
FOLHETO ILLUSTRATIVO
ΦΥΛΛΑΔΙΟ ΕΠΕΞΗΓΗΣΕΩΝ
説明書

* Refers to Class IIa products only
HPF S.p.A. CE Certification

Date: March 23, 2011

English	5	EN
Français	13	FR
Deutsch	21	DE
Dansk	29	DA
Nederlands	37	NL
Español	45	ES
Italiano	53	IT
Svenska	61	SE
Português	69	PT
Ελληνικά	77	GR
日本語	85	JP

ENGLISH
INFORMATION LEAFLET
HPF S.p.A. INSTRUMENTS

All. 9 B_1 Ed.1 del 27-01-12

WARNINGS

- Critical elements, such as handles, articulated joint instruments or threads and in particular all those medical instruments with atraumatic toothings require particularly careful inspections.
- In order to avoid damage and subsequent corrosion due to rubbing of metals, avoid under all circumstances the use of metal brushes or metal sponges for removing stains.
- Do not exceed 137 °C.
- Devices that cannot be disassembled or that are complex require particular attention during cleaning.
- After having disinfected and cleaned the instruments, immediately dry them.
- In case of thermolabile material, steam sterilization is preferable.
- Machine treatment with thermal disinfection and steam sterilization should be applied.
- Acetabular cutters: do not use hand-pieces with a torque higher than 30 Nm and the axial thrust is allowed to a single user.
- Flexible mandrel: do not subject the instrument to a torque higher than 7 Nm.

DESCRIPTION

HPF S.p.A. instruments are medical instruments used for surgeries:

- manual medical instruments belonging to class I;
- energy transfer instruments belonging to class II a (Refer to the Directive 93/42/EEC).

HPF S.p.A. medical devices are authorized for use only by specialized, trained personnel certified to perform surgeries or handle related surgical instruments.

USE

Before clinical use, the surgeon must be familiar with all aspects of the surgery and the limits of the instruments.

These instruments have a duration time limit; this limit is determined by the level of wear and tear and damage arising from the intended use.

Therefore, it is recommended to replace the instruments when the user notices reduced performance or in the case visible signs of wear and tear appear.

RECOMMENDATIONS FOR USE

Generally, device treatment includes:

- Preparation (pre-treatment, gathering of instruments, pre-cleaning and disassembly, if applicable).
- Cleaning, disinfection, rinsing and drying.
- Visual inspection of cleaning and perfect conditions of the material.
- Verification of correct operation.

- Packaging and sterilization, authorization and storing, whenever applicable.

After use, the instruments must be properly conditioned.

MANUAL CLEANING

- Reusable instruments must be cleaned immediately after use.
- Recesses (cavities) and hidden areas within an instrument must be inspected periodically in order to ensure the absence of residual substances.
- For manual cleaning, we recommend using a soft cloth with low particle release or synthetic brushes.
- After manual disinfectant cleaning it is necessary to perform intense and thorough rinsing using clear running water. During the procedure, manually remove any possible dirt residue. Then proceed to a final rinse using demineralized water and thoroughly dry the instruments. A good cleaning of the instruments plays also an important role in maintaining the value of the instruments and sets the stage for a successful sterilization.
- Pursuant to the international standards (EN ISO 15883) and national guidelines, only validated procedures for cleaning and disinfecting should be applied.

MECHANICAL CLEANING

- Reusable instruments must be cleaned immediately after use.
- Recesses (cavities) and hidden areas within an instrument must be inspected periodically in order to ensure the absence of residual substances.
- The best method to achieve standardization, cleaning and disinfection is to use the machine treatment method. A good cleaning of the instruments plays also an important role in maintaining the value of the instruments and sets the stage for a successful sterilization.
- Pursuant to the international standards (EN ISO 15883) and national guidelines, only validated procedures for cleaning and disinfecting should be applied.

STERILIZATION

- This instrument exits the premises of **HPF S.p.A.** non-sterilized.
- All instruments and parts of instruments must be sterilized before use.
- The use of sterilized instruments on the patient establishes, within the framework of the European regulations in force that all the instruments are adequately cleaned and disinfected, sterilized in a suitable sterile packaging using a validated procedure and stored according to the regulations in

force for sterile materials. This is why it is important to use only methods and sterilization devices that allow a validated sterilization process.

- All instruments must be removed from their packaging prior to sterilization.
- Reusable instruments must be disinfected and sterilized after use.
- The recommended sterilization for the instruments manufactured by **HPF** S.p.A. is steam sterilization (autoclave) according to the EN 285 Directive and European pharmacopoeia regulations:
 - **HPF** S.p.A. accepts no responsibility in case problems arising from sterilizations different from the recommended one.
 - The instruments must not come into contact with fluoridated or chloric treated products or with greasy detergents.
 - Instruments with synthetic material elements (plastics) must never come in contact with strong acids (pH'4), with bases, organic solvents or ammonia, and not even with oxidized chemical products or any type product that may alter the materials.

PROHIBITION

It is strictly forbidden to apply any modification to the instruments manufactured by **HPF** S.p.A.;

make such modifications. In case of failure to heed this prohibition, **HPF** S.p.A. accepts no responsibility for the consequences.

STORAGE

- The instruments are sensitive to scratches and denting. Scratches on their surface can increase the risk of corrosion. For this reason it is recommended to handle them with care.
- Storage of medical instruments must be in an environment free from moisture in order to avoid excessive corrosion. This is valid also for the transportation and conditioning of the instruments.
- The surgical instruments are subject to wear and tear; thus, prior to use it is necessary to inspect them to verify that they are free from defects.

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